Oxoid™ Prepared Medium

Oxoid Bordetella Selective Medium

REF PB5065A

Intended Use

IVD

Oxoid Bordetella Selective Medium is a selective medium for the cultivation and isolation of fastidious organisms especially *Bordetella pertussis* and *Bordetella parapertussis* from nasopharyngeal/pernasal samples.

The device is for professional use only, is not automated, and nor is it a companion diagnostic.

Summary and Explanation

Bordetella spp. are Gram-negative, aerobic coccobacilli which mainly inhabit the surface of the respiratory tract of man and other warm-blooded animals. *B. pertussis* usually infects and causes severe respiratory disease (commonly known as whooping cough) in young children, although it may cause disease in adults. *B. parapertussis* is also responsible for respiratory tract infections in man while these tend to be milder, they may be severe in immunocompromised patients.

Principle of Method

Charcoal Agar used in the preparation of Oxoid Bordetella Selective Medium was developed by Oxoid for the cultivation of *B. pertussis* and *Haemophilus influenzae*. Proom¹ showed that nicotinic acid was an essential growth factor for the bordetellae. Ensminger *et al.*² used a charcoal medium for the growth of *Bordetella pertussis* in vaccine production and found that the medium could replace Bordet-Gengou. Mishulow *et al.*³ used charcoal agar for *B. pertussis* cultivation. Bradford and Slavin⁴ concluded that poor results with the cultivation of post-nasal swabs for *B. pertussis* were due to overgrowth by competing microorganisms rather than the absence of *B. pertussis*. Lacey⁵ was able to confirm this. He noted that even with the addition of 0.25 IU/mL of penicillin⁶, the detection of *B. pertussis* was often made tedious and difficult because of overgrowth by penicillin insensitive flora. The fact that this occurred much more often with post-nasal rather than per-nasal swabs was felt to account for the relatively small use made of per-nasal swabs in the diagnosis of whooping cough.

Lacey⁵ then found that addition of M&B 938 (4,4'-diamidino-diphenylamine dihydrochloride), at a concentration of 12 µg/mL in a partially defined agar base, together with penicillin, increased the selectivity of the medium by partially inhibiting the growth of penicillin insensitive strains of *Haemophilus*, *Staphylococcus* and *Neisseria* and led to a higher isolation rate for *B. pertussis*.

Addition of 40 µg/mL cephalaxin to the medium⁷ instead of penicillin and M&B 938 makes it more selective and significantly increases the isolation rate of *Bordetella* species. Coliforms are inhibited and many microorganisms are totally suppressed; some exceptions are *Pseudomonas aeruginosa* and fungi; for suppression of fungi the addition of 50 µg/mL of amphotericin B is recommended. The use of cephalaxin instead of penicillin allows a greater recovery of stressed *B. pertussis* cells thus further increasing the isolation rate.

Typical Formula

	grams per litre
'Lab-Lemco' powder	10.0
Peptone	10.0
Starch	10.0
Bacteriological charcoal	4.0
Sodium chloride	5.0
Nicotinic acid	0.001
Agar	12.0

Additions

	per litre
Cephalaxin	40 mg
Defibrinated Horse Blood	100.0 ml

Physical Characteristics

Colour	Jet black
Clarity	Opaque
Fill weight	17 g ± 5 %
pH	7.5 ± 0.2

Materials Provided

- The pack contains 10 x 90mm agar plates, film wrapped.
- Each plate should only be used once.
- Each pack contains enough plates for 10 individual tests.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms

Storage

This product is ready to use and no further preparation is necessary. Store product in its original packaging at 2–12°C until used. Store away from light. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

Warnings and Precautions

This product is for *in vitro* diagnostic use and should only be used by trained individuals. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully.

Inspect the product packaging before first use. Do not use the product if there is any visible damage to the packaging or plates. Do not use the product beyond the stated expiry date. Do not use the device if signs of contamination are present. Do not use the device if the colour has changed or there are other signs of deterioration.

For professional use only. Refer to the Safety Data Sheet (SDS) at www.thermofisher.com

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following the recommended guidelines.

Procedure

Refer to local protocols and guidelines.

- Inoculate the specimen onto the medium and streak to obtain isolated colonies.
- If sample material is being cultured directly from a swab, roll the swab over a small area and streak to obtain single colonies.
- Incubate the plates at 35°C for 3 days in a moist chamber. Do not discard the plates as negative until they have been incubated for 6 days.

B. pertussis colonies are almost transparent with a 'bisected' pearl appearance.

Quality Control

This medium can be tested with the following strains:

Incubation Conditions: 5 days @ 36 ± 1°C, aerobic, humid atmosphere.

Positive Controls	
Inoculum of 50 – 120 colony forming units (cfu).	
<i>Bordetella parapertussis</i> ATCC®15311	1 – 2 mm, grey shiny colonies
Colony count shall be greater than or equal to 50% of the control medium (Bordetella selective medium).	
Inoculum of 10 ³ – 10 ⁴ cfu	
<i>Bordetella pertussis</i> ATCC®9797	1 mm, grey shiny colonies
Negative Control	
Inoculum of ≥ 10 ⁴ cfu	
<i>Staphylococcus aureus</i> ATCC®25923	Complete inhibition (≤ 10 cfu)

Note:

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

Performance

Performance was evaluated at 36 ± 1°C in an aerobic, humid atmosphere using 7 bacterial strains *B. pertussis* (3) and one each of *B. parapertussis*, *Haemophilus influenzae*, *Staphylococcus aureus* and *Klebsiella pneumoniae*. Strains of *B. pertussis* and *B. parapertussis* gave good growth within 5 days. *Haemophilus influenzae*, *S. aureus* and *K. pneumoniae* were inhibited.

Limitations

Coliforms are inhibited and many microorganisms are totally suppressed however *Pseudomonas aeruginosa* and fungi may grow through. Some strains of *Bordetella* spp. with particular nutritional requirements may not grow on this medium.

Bibliography

1. Proom H. (1955) *J. Gen. Microbiol.* 12 (1). 63-75.
2. Ensminger P.W., Culbertson C.G. and Powell H.M. (1953) *J. Infect. Dis.* 93 (3). 266-268.
3. Mishulow Lucy, Sharpe L.S. and Choen Lillian L. (1953) *Amer. J. Pub. Health* 43 (11). 1466-1472
4. Bradford W. L. and Slavin B. (1940) *Proc. Soc. Exp. Biol.* 43. 590-3.
5. Lacey B. W. (1954) *J. Hyg.* 52 (3) 273 - 303.
6. Fleming A. (1932) *J. Path. Bact.* 35. 831 - 842
7. Regan J. and Lowe F. (1977) *J. Clin. Microbiol.* 6. 303-309

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device

	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Manufacturer
	Authorized representative in the European Community/European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Unique device identifier

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. ATCC® is a trademark of American Type Culture Collection. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

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For technical assistance please contact your local distributor.

Version	Date of issue and modifications introduced
1.0	2022-08-04. New document

Oxoid™ Prepared Medium

Oxoid Bordetella Selective Medium

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Oxoid Deutschland GmbH, Am Lippeglacis 4-8, 46483 Wesel, Alemania